

K001069 MODEL 4000 EEG MONITOR WITH PSI (PSA4000)Jun 30, 2000
88 days to decisionK001069 · Product code: **OLW** · Neurology
Source: <https://www.510kdatabase.net/k001069/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Index-generating Electroencephalograph Software (OLW)
Date received	Apr 3, 2000
Decision date	Jun 30, 2000
Days to decision	88 days
Third-party review	No
Summary / Statement	Summary

APPLICANT

Company	Physiometrix, Inc.
Location	Sunnyvale, CA, US
Contact	DAWN E FRAZER
510(k) history	22 submissions · 22 cleared · 1992-2005

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k001069/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated June 29, 2026