

**K002257 TECNO DISPOSABLE FINGERSWITCH PENCIL,
MODEL 150-100**Oct 19, 2000
86 days to decisionK002257 · Product code: **GEI** · General & Plastic SurgerySource: <https://www.510kdatabase.net/k002257/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Electrosurgical, Cutting & Coagulation & Accessories (GEI)
Date received	Jul 25, 2000
Decision date	Oct 19, 2000
Days to decision	86 days
Third-party review	No
Summary / Statement	Summary

APPLICANT

Company	Tecno Instruments (Pvt.) , Ltd.
Location	Great Neck, NY, US
Contact	SUSAN D GOLDSTEIN-FALK
510(k) history	2 submissions · 2 cleared · 2000-2004

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k002257/>; Data sourced from FDA 510(k) public records (accessdata.fda.gov).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated July 3, 2026