

**K003444 MODEL WS-500 DIGITAL WRIST BLOOD PRESSURE MONITOR**

Nov 21, 2000  
15 days to decision

K003444 · Product code: **DXN** · Cardiovascular  
Source: <https://www.510kdatabase.net/k003444/>

**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	System, Measurement, Blood-pressure, Non-invasive (DXN)
Date received	Nov 6, 2000
Decision date	Nov 21, 2000
Days to decision	15 days
Third-party review	Yes
Summary / Statement	Summary

**APPLICANT**

---

Company	<b>Nihon Seimitsu Sokki Co., Ltd.</b>
Location	Shibukawa-Shi, Gunma-Ken, JP
Contact	ROBERT D WAXHAM
510(k) history	21 submissions · 21 cleared · 1985-2012

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)

Device record: <https://www.510kdatabase.net/k003444/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated June 14, 2026