

**K022281 SAFE MAXI VENOUS/CARDIOTOMY RESERVIOR,
MODEL 016045**Aug 6, 2002
22 days to decisionK022281 · Product code: **DTN** · Cardiovascular
Source: <https://www.510kdatabase.net/k022281/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Reservoir, Blood, Cardiopulmonary Bypass (DTN)
Date received	Jul 15, 2002
Decision date	Aug 6, 2002
Days to decision	22 days
Third-party review	No
Summary / Statement	Summary

APPLICANT

Company	Polystan A/S
Location	Denmark, DK
Contact	DANA OLSEN
510(k) history	6 submissions · 6 cleared · 1991-2002

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k022281/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated June 28, 2026