

**K030038 PHILIPS MP60, MP70 AND MP90 INTELLIVUE PATIENT MONITORS WITH PORTAL TECHNOLOGY AND M3185 OPT C41 WIRELESS BEDSIDE ADAPTR**Mar 12, 2003  
65 days to decisionK030038 · Product code: **MHX** · Cardiovascular  
Source: <https://www.510kdatabase.net/k030038/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Special
Device classification	Monitor, Physiological, Patient(with Arrhythmia Detection Or Alarms) (MHX)
Date received	Jan 6, 2003
Decision date	Mar 12, 2003
Days to decision	65 days
Third-party review	No
Summary / Statement	Summary

**APPLICANT**

---

Company	<b>Philips Medical Systems, Inc.</b>
Location	Mchenry, IL, US
Contact	DAVID OSBORN
510(k) history	111 submissions · 110 cleared · 1983-2010

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k030038/>; Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated May 17, 2026