

**K030077 GISH TUBING AND CONNECTORS WITH GBS  
COATING**Mar 14, 2003  
65 days to decisionK030077 · Product code: **DTL** · Cardiovascular  
Source: <https://www.510kdatabase.net/k030077/>**SUBMISSION DETAILS**

---

|                       |  |
|-----------------------|--|
| Decision              | Substantially Equivalent (Cleared)                                 |
| Submission type       | Traditional  |
| Device classification | Adaptor, Stopcock, Manifold, Fitting, Cardiopulmonary Bypass (DTL) |
| Date received         | Jan 8, 2003  |
| Decision date         | Mar 14, 2003   |
| Days to decision      | 65 days  |
| Third-party review    | No   |
| Summary / Statement   | Summary  |

**APPLICANT**

---

|                |   |
|----------------|---|
| Company        | <b>Gish Biomedical, Inc.</b>            |
| Location       | Mchenry, IL, US                         |
| Contact        | J. HARVEY KNAUSS                        |
| 510(k) history | 75 submissions · 75 cleared · 1983-2009 |

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k030077/>. Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 26, 2026