

**K040304 PHILIPS INTELLIVUE PATIENT MONITOR, MODELS  
MP40 & MP50**Mar 1, 2004  
21 days to decisionK040304 · Product code: **MHX** · Cardiovascular  
Source: <https://www.510kdatabase.net/k040304/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Special
Device classification	Monitor, Physiological, Patient(with Arrhythmia Detection Or Alarms) (MHX)
Date received	Feb 9, 2004
Decision date	Mar 1, 2004
Days to decision	21 days
Third-party review	No
Summary / Statement	Summary

**APPLICANT**

---

Company	<b>Philips Medizin Systeme Boeblingen GmbH</b>
Location	B?blingen, DE
Contact	HAUKE SCHIK
510(k) history	48 submissions · 48 cleared · 2004-2026

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k040304/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 18, 2026