

K041235 PHILIPS MP20, MP30, MP40, MP50. MP60, MP70 AND MP90 INTELLIVUE PATIENT MONITORSJun 4, 2004
24 days to decisionK041235 · Product code: **MHX** · Cardiovascular
Source: <https://www.510kdatabase.net/k041235/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Special
Device classification	Monitor, Physiological, Patient(with Arrhythmia Detection Or Alarms) (MHX)
Date received	May 11, 2004
Decision date	Jun 4, 2004
Days to decision	24 days
Third-party review	No
Summary / Statement	Summary

APPLICANT

Company	Philips Medizin Systeme Boeblingen GmbH
Location	B?blingen, DE
Contact	MARKUS STACHA
510(k) history	48 submissions · 48 cleared · 2004-2026

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k041235/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 18, 2026