

**K050141 THE PHILIPS MP20, MP30, MP40, MP50, MP60, MP70  
AND MP90 INTELLIVUE PATIENT MONITORS**Feb 11, 2005  
17 days to decisionK050141 · Product code: **MHX** · Cardiovascular  
Source: <https://www.510kdatabase.net/k050141/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Special
Device classification	Monitor, Physiological, Patient(with Arrhythmia Detection Or Alarms) (MHX)
Date received	Jan 25, 2005
Decision date	Feb 11, 2005
Days to decision	17 days
Third-party review	No
Summary / Statement	Summary

**APPLICANT**

---

Company	<b>Philips Medizin Systeme Boeblingen GmbH</b>
Location	B?blingen, DE
Contact	HAUKE SCHIK
510(k) history	48 submissions · 48 cleared · 2004-2026

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k050141/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 18, 2026