

**K060963 TWIN STAR COMPARTMENT PRESSURE MONITORING AND FLUID COLLECTION MONITOR (CMS MONITOR)**Jun 6, 2006  
60 days to decisionK060963 · Product code: **LXC** · Orthopedic  
Source: <https://www.510kdatabase.net/k060963/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Monitor, Pressure, Intracompartmental (LXC)
Date received	Apr 7, 2006
Decision date	Jun 6, 2006
Days to decision	60 days
Third-party review	No
Summary / Statement	Summary

**APPLICANT**

---

Company	<b>Twin Star Medical</b>
Location	St. Paul, MN, US
Contact	CHET SIEVERT
510(k) history	1 submissions · 1 cleared · 2006-2006

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k060963/> Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated May 19, 2026