

**K073445 CLEVER TD-3026 DIGITAL WRIST BLOOD PRESSURE MONITOR, MODEL TD-3026**Apr 2, 2008  
117 days to decisionK073445 · Product code: **DXN** · Cardiovascular  
Source: <https://www.510kdatabase.net/k073445/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Special
Device classification	System, Measurement, Blood-pressure, Non-invasive (DXN)
Date received	Dec 7, 2007
Decision date	Apr 2, 2008
Days to decision	117 days
Third-party review	No
Summary / Statement	Statement

**APPLICANT**

---

Company	<b>Taidoc Technology Corporation</b>
Location	Taipei, County, TW
Contact	ERICA LI
510(k) history	123 submissions · 123 cleared · 2004-2024

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k073445/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated June 25, 2026