

K081590 YA HORNG ARM TYPE DIGITAL BLOOD PRESSURE MONITOR, MODELS BP-600, BP-600R, BP-600U, BP-600B, BP-600RB, BP-600UB, BP-600JJun 19, 2008
13 days to decisionK081590 · Product code: **DXN** · Cardiovascular
Source: <https://www.510kdatabase.net/k081590/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	System, Measurement, Blood-pressure, Non-invasive (DXN)
Date received	Jun 6, 2008
Decision date	Jun 19, 2008
Days to decision	13 days
Third-party review	Yes
Summary / Statement	Summary

APPLICANT

Company	Ya Horng Electronic Co., Ltd.
Location	Hsin-Chu City, TW
Contact	STEVEN CHANG
510(k) history	10 submissions · 10 cleared · 2004-2022

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k081590/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 26, 2026