

K092045 FULLY AUTOMATIC ELECTRONIC BLOOD PRESSURE MONITOR, MODELS KD-738 AND KD-739Jul 30, 2009
24 days to decisionK092045 · Product code: **DXN** · Cardiovascular
Source: <https://www.510kdatabase.net/k092045/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Special
Device classification	System, Measurement, Blood-pressure, Non-invasive (DXN)
Date received	Jul 6, 2009
Decision date	Jul 30, 2009
Days to decision	24 days
Third-party review	No
Summary / Statement	Summary

APPLICANT

Company	Andon Health Co, Ltd.
Location	Tiajin, CN
Contact	LIU YI
510(k) history	92 submissions · 92 cleared · 2008-2025

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k092045/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated June 22, 2026