

**K092179 APLIO XG DIAGNOSTIC ULTRASOUND SYSTEM,  
MODEL SSA - 790A VERSION 4.0R001**Jul 29, 2009  
8 days to decisionK092179 · Product code: IYN · Radiology  
Source: <https://www.510kdatabase.net/k092179/>**SUBMISSION DETAILS**

---

|                       |   |
|-----------------------|---|
| Decision              | Substantially Equivalent (Cleared)                |
| Submission type       | Traditional                                       |
| Device classification | System, Imaging, Pulsed Doppler, Ultrasonic (IYN) |
| Date received         | Jul 21, 2009                                      |
| Decision date         | Jul 29, 2009                                      |
| Days to decision      | 8 days  |
| Third-party review    | Yes   |
| Summary / Statement   | Summary   |

**APPLICANT**

---

|                |  |
|----------------|--|
| Company        | <b>Toshiba America Medical Systems, In.C</b> |
| Location       | Tustin, CA, US                               |
| 510(k) history | 146 submissions · 146 cleared · 1989-2014    |

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k092179/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 16, 2026