

**K100207 NON-MYDRIATIC RETINAL CAMERA, MODEL TRC-NW8F**Feb 12, 2010  
18 days to decisionK100207 · Product code: **HKI** · Ophthalmic  
Source: <https://www.510kdatabase.net/k100207/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Camera, Ophthalmic, Ac-powered (HKI)
Date received	Jan 25, 2010
Decision date	Feb 12, 2010
Days to decision	18 days
Third-party review	Yes
Summary / Statement	Summary

**APPLICANT**

---

Company	<b>Topcon Corp.</b>
Location	Paramus, NJ, US
Contact	TAKAO SUGAWARA
Website	<a href="http://www.topcon.com">http://www.topcon.com</a>
510(k) history	9 submissions · 9 cleared · 2004-2014

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k100207/> Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated June 6, 2026