

**K101681 TRANSTEK BLOOD PRESSURE MONITOR MODEL  
TMB-986, TMB-987, TMB-995**May 27, 2011  
346 days to decisionK101681 · Product code: **DXN** · Cardiovascular  
Source: <https://www.510kdatabase.net/k101681/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Abbreviated
Device classification	System, Measurement, Blood-pressure, Non-invasive (DXN)
Date received	Jun 15, 2010
Decision date	May 27, 2011
Days to decision	346 days
Third-party review	No
Summary / Statement	Summary

**APPLICANT**

---

Company	<b>Zhongshan Transtek Electronics Co., Ltd.</b>
Location	Chengdu, Sichuan, CN
Contact	LEO WANG
510(k) history	9 submissions · 9 cleared · 2011-2013

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k101681/> Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated June 9, 2026