

K102988 ELECTRONIC BLOOD PRESURE MONITORNov 9, 2010
34 days to decisionK102988 · Product code: **DXN** · Cardiovascular
Source: <https://www.510kdatabase.net/k102988/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	System, Measurement, Blood-pressure, Non-invasive (DXN)
Date received	Oct 6, 2010
Decision date	Nov 9, 2010
Days to decision	34 days
Third-party review	Yes
Summary / Statement	Summary

APPLICANT

Company	Shenzhen Pango Electronic Co., Ltd.
Location	Shanghai, CN
Contact	Diana Hong
510(k) history	9 submissions · 9 cleared · 2010-2021

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k102988/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 26, 2026