

**K122106 SPIN DRIVE ALWAYS-ON TIP TRACKED STEERABLE  
CATHETER ALWAYS-ON TIP TRACKED SHEATH ALWAYS-ON  
TIP TRACKED ASPIRATION NEEDLE**Dec 21, 2012  
157 days to decisionK122106 · Product code: **JAK** · Radiology  
Source: <https://www.510kdatabase.net/k122106/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	System, X-ray, Tomography, Computed (JAK)
Date received	Jul 17, 2012
Decision date	Dec 21, 2012
Days to decision	157 days
Third-party review	No
Summary / Statement	Summary

**APPLICANT**

---

Company	<b>Veran Medical Technologies, Inc.</b>
Location	Nashville, TN, US
Contact	SCOTT WRIGHTSTONE
510(k) history	7 submissions · 7 cleared · 2006-2021

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k122106/> Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated May 27, 2026