

**K142951 MD210C3 21.3 Diagnostic Imaging LCD Monitor**Dec 10, 2014  
61 days to decisionK142951 · Product code: **PGY** · Radiology  
Source: <https://www.510kdatabase.net/k142951/>**SUBMISSION DETAILS**

---

|                       |                                     |
|-----------------------|-------------------------------------|
| Decision              | Substantially Equivalent (Cleared)  |
| Submission type       | Traditional                         |
| Device classification | Display, Diagnostic Radiology (PGY) |
| Date received         | Oct 10, 2014                        |
| Decision date         | Dec 10, 2014                        |
| Days to decision      | 61 days                             |
| Third-party review    | No                                  |
| Summary / Statement   | Summary                             |

**APPLICANT**

---

|                |   |
|----------------|---|
| Company        | <b>Nec Display Solutions , Ltd.</b>     |
| Location       | Ashigarakami-Gun, Kanagawa, JP          |
| Contact        | Satoru Kotani                           |
| 510(k) history | 15 submissions · 15 cleared · 2006-2016 |

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k142951/> Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated May 16, 2026