

K150430 Arm Blood Pressure Monitor, Wrist Blood Pressure MonitorApr 17, 2015
57 days to decisionK150430 · Product code: **DXN** · Cardiovascular
Source: <https://www.510kdatabase.net/k150430/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	System, Measurement, Blood-pressure, Non-invasive (DXN)
Date received	Feb 19, 2015
Decision date	Apr 17, 2015
Days to decision	57 days
Third-party review	No
Summary / Statement	Summary

APPLICANT

Company	Fudakang Industrial Co., Ltd.
Location	Buena Park, CA, US
Contact	BOB YU
510(k) history	4 submissions · 4 cleared · 2012-2018

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k150430/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 30, 2026