

**K160347 Hisense LDC monitor models
HMD2G21/HMD3G21/HMD5G21**Mar 3, 2016
24 days to decisionK160347 · Product code: **PGY** · Radiology
Source: <https://www.510kdatabase.net/k160347/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Display, Diagnostic Radiology (PGY)
Date received	Feb 8, 2016
Decision date	Mar 3, 2016
Days to decision	24 days
Third-party review	No
Summary / Statement	Summary

APPLICANT

Company	Qingdao Hisense Medical Equipment Co., Ltd.
Location	Qingdao, CN
Contact	Liu zhitang
510(k) history	9 submissions · 9 cleared · 2015-2024

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k160347/>. Data sourced from FDA 510(k) public records (accessdata.fda.gov).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 26, 2026