

K161712 Arm automatic blood pressure monitorApr 5, 2017
288 days to decisionK161712 · Product code: **DXN** · Cardiovascular
Source: <https://www.510kdatabase.net/k161712/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Special
Device classification	System, Measurement, Blood-pressure, Non-invasive (DXN)
Date received	Jun 21, 2016
Decision date	Apr 5, 2017
Days to decision	288 days
Third-party review	No
Summary / Statement	Summary

APPLICANT

Company	Shenzhen Pump Medical System Co., Ltd.
Location	Shenzhen, Guangdong, CN
Contact	XIE QIONGYU
510(k) history	7 submissions · 7 cleared · 2013-2017

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k161712/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated June 1, 2026