

**K161794 CENTURION Vision System (Active Sentry™)**Apr 14, 2017  
288 days to decisionK161794 · Product code: **HQC** · Ophthalmic  
Source: <https://www.510kdatabase.net/k161794/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Unit, Phacofragmentation (HQC)
Date received	Jun 30, 2016
Decision date	Apr 14, 2017
Days to decision	288 days
Third-party review	No
Summary / Statement	Summary

**APPLICANT**

---

Company	<b>Alcon Research, Ltd.</b>
Location	Fort Worth, TX, US
Contact	Paul Swift
510(k) history	16 submissions · 16 cleared · 2000-2017

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k161794/> Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated June 1, 2026