

**K170527 Phoenix Clinical ICON**Sep 15, 2017  
205 days to decisionK170527 · Product code: **HKI** · Ophthalmic  
Source: <https://www.510kdatabase.net/k170527/>**SUBMISSION DETAILS**

---

|                       |                                      |
|-----------------------|--------------------------------------|
| Decision              | Substantially Equivalent (Cleared)   |
| Submission type       | Traditional                          |
| Device classification | Camera, Ophthalmic, Ac-powered (HKI) |
| Date received         | Feb 22, 2017                         |
| Decision date         | Sep 15, 2017                         |
| Days to decision      | 205 days                             |
| Third-party review    | No                                   |
| Summary / Statement   | Summary                              |

**APPLICANT**

---

|                |   |
|----------------|---|
| Company        | <b>Phoenix Technology Group, Inc. Db a Phoenix Clinical</b> |
| Location       | Pleasanton, CA, US  |
| Contact        | Christopher A. Henderson                                    |
| 510(k) history | 1 submissions · 1 cleared · 2017-2017                       |

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k170527/> Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated May 31, 2026