

**K182891 TEEClean(R) Automated TEE Probe Cleaner  
Disinfectant with TEEZyme(R) Cleaner and TD-5(R) or TD-8(R)  
High Level Disinfectants**Apr 12, 2019  
179 days to decisionK182891 · Product code: **PSW** · General Hospital  
Source: <https://www.510kdatabase.net/k182891/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	High Level Disinfection Reprocessing Instrument For Ultrasonic Transducers, Liquid (PSW)
Date received	Oct 15, 2018
Decision date	Apr 12, 2019
Days to decision	179 days
Third-party review	No
Summary / Statement	Summary

**APPLICANT**

---

Company	<b>Cs Medical, LLC</b>
Location	Raleigh, NC, US
Contact	Kendall Ashe
510(k) history	6 submissions · 6 cleared · 2005-2023

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k182891/> Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated May 26, 2026