

K791233 MODEL EEG-4214 ELECTROENCEPHALOGRAPHSep 12, 1979
75 days to decisionK791233 · Product code: **GWQ** · Neurology
Source: <https://www.510kdatabase.net/k791233/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Full-montage Standard Electroencephalograph (GWQ)
Date received	Jun 29, 1979
Decision date	Sep 12, 1979
Days to decision	75 days
Third-party review	No

APPLICANT

Company	Nihon Kohden America, Inc.
Location	Foothill Ranch, CA, US
510(k) history	166 submissions · 163 cleared · 1979-2012

510k Database - www.510kdatabase.net

Device record: <https://www.510kdatabase.net/k791233/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 24, 2026