

**K801064 DISPOSABLE PATIENT CABLE MODEL 5807**May 14, 1980  
8 days to decisionK801064 · Product code: **DSA** · CardiovascularSource: <https://www.510kdatabase.net/k801064/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Cable, Transducer And Electrode, Patient, (including Connector) (DSA)
Date received	May 6, 1980
Decision date	May 14, 1980
Days to decision	8 days
Third-party review	No

**APPLICANT**

---

Company	<b>Daig Corp.</b>
Location	Mchenry, IL, US
510(k) history	63 submissions · 63 cleared · 1977-2000

---

**510k Database** - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k801064/> Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated July 3, 2026