

**K851400 SIRECUST 420 PATIENT DIALOG MONITOR**Jun 26, 1985  
78 days to decisionK851400 · Product code: **DXJ** · Cardiovascular  
Source: <https://www.510kdatabase.net/k851400/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Display, Cathode-ray Tube, Medical (DXJ)
Date received	Apr 9, 1985
Decision date	Jun 26, 1985
Days to decision	78 days
Third-party review	No

**APPLICANT**

---

Company	<b>Siemens Medical Solutions USA, Inc.</b>
Location	Hoffman Estates, IL, US
Contact	HARRY K SCHWILL
510(k) history	778 submissions · 778 cleared · 1980-2026

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k851400/>; Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated May 15, 2026