

**K860531 SONOLAYER-L, SAL-38A, WITH SCANNER SAL-38AS**Jul 18, 1986  
156 days to decisionK860531 · Product code: **ITX** · Radiology  
Source: <https://www.510kdatabase.net/k860531/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Transducer, Ultrasonic, Diagnostic (ITX)
Date received	Feb 12, 1986
Decision date	Jul 18, 1986
Days to decision	156 days
Third-party review	No

**APPLICANT**

---

Company	<b>Toshiba Medical Systems</b>
Location	Mchenry, IL, US
Contact	WYATT CANNADY
510(k) history	54 submissions · 54 cleared · 1977-1996

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k860531/>, Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated May 20, 2026