

**K863243 SHINER STRESS DIRECTING ATTACHMENT**Oct 30, 1986  
69 days to decisionK863243 · Product code: **EGG** · Dental  
Source: <https://www.510kdatabase.net/k863243/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Attachment, Precision, All (EGG)
Date received	Aug 22, 1986
Decision date	Oct 30, 1986
Days to decision	69 days
Third-party review	No

**APPLICANT**

---

Company	<b>Js &amp; R - R &amp; D Co.</b>
Location	San Diego, CA, US
Contact	JIM SHINER
510(k) history	1 submissions · 1 cleared · 1986-1986

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k863243/>; Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated July 5, 2026