

**K870825 NEW MODEL STYLES FOR AORTIC ARCH CANNULA**Mar 23, 1987  
20 days to decisionK870825 · Product code: **DQR** · CardiovascularSource: <https://www.510kdatabase.net/k870825/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Cannula, Catheter (DQR)
Date received	Mar 3, 1987
Decision date	Mar 23, 1987
Days to decision	20 days
Third-party review	No

**APPLICANT**

---

Company	<b>Shiley, Inc.</b>
Location	Mchenry, IL, US
Contact	MORTON BARAK
510(k) history	174 submissions · 174 cleared · 1976-1993

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k870825/> Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated May 26, 2026