

K872517 TL 400 TOTALAB VASCULAR DIAGNOSTIC SYSTEMOct 1, 1987
99 days to decisionK872517 · Product code: **JOM** · CardiovascularSource: <https://www.510kdatabase.net/k872517/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Plethysmograph, Photoelectric, Pneumatic Or Hydraulic (JOM)
Date received	Jun 24, 1987
Decision date	Oct 1, 1987
Days to decision	99 days
Third-party review	No

APPLICANT

Company	D. E. Hokanson, Inc.
Location	Mchenry, IL, US
Contact	EUGENE HOKANSON
510(k) history	19 submissions · 19 cleared · 1979-2004

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k872517/>; Data sourced from FDA 510(k) public records (accessdata.fda.gov).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 19, 2026