

**K890687 ECT-100 ELECTRONIC CORNEAL TOPOGRAPHER**Sep 5, 1989  
208 days to decisionK890687 · Product code: **HLQ** · Ophthalmic  
Source: <https://www.510kdatabase.net/k890687/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Keratoscope, Ac-powered (HLQ)
Date received	Feb 9, 1989
Decision date	Sep 5, 1989
Days to decision	208 days
Third-party review	No

**APPLICANT**

---

Company	<b>Optimed, Inc.</b>
Location	Oakland, CA, US
Contact	DOUGLAS DONALDSON
510(k) history	4 submissions · 4 cleared · 1988-1990

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k890687/>; Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated June 30, 2026