

K897063 FUKUDA DENSHI MODEL DS-3300 ETRFeb 1, 1990
43 days to decisionK897063 · Product code: **DSI** · Cardiovascular
Source: <https://www.510kdatabase.net/k897063/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Detector And Alarm, Arrhythmia (DSI)
Date received	Dec 20, 1989
Decision date	Feb 1, 1990
Days to decision	43 days
Third-party review	No

APPLICANT

Company	Fukuda Denshi USA, Inc.
Location	Mchenry, IL, US
Contact	ROBERT J STEURER
510(k) history	68 submissions · 68 cleared · 1984-2018

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k897063/>; Data sourced from FDA 510(k) public records (accessdata.fda.gov).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 18, 2026