

**K913193 MODEL 100 CENTRAL STATION**Oct 16, 1991  
90 days to decisionK913193 · Product code: **DQK** · CardiovascularSource: <https://www.510kdatabase.net/k913193/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Computer, Diagnostic, Programmable (DQK)
Date received	Jul 18, 1991
Decision date	Oct 16, 1991
Days to decision	90 days
Third-party review	No
Summary / Statement	Statement

**APPLICANT**

---

Company	<b>Protocol Systems, Inc.</b>
Location	Beaverton, OR, US
Contact	JAMES W SANDBERT
510(k) history	15 submissions · 14 cleared · 1988-2000

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k913193/> Data sourced from FDA 510(k) public records ([accessdata.fda.gov](http://accessdata.fda.gov)).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at [accessdata.fda.gov](http://accessdata.fda.gov). - Generated May 26, 2026