

K935103 VISIONLINK ACR-NEMA CONVERSION UTILITIESFeb 9, 1994
107 days to decisionK935103 · Product code: **LMD** · Radiology
Source: <https://www.510kdatabase.net/k935103/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	System, Digital Image Communications, Radiological (LMD)
Date received	Oct 25, 1993
Decision date	Feb 9, 1994
Days to decision	107 days
Third-party review	No
Summary / Statement	Summary

APPLICANT

Company	Evergreen Technologies, Inc.
Location	Gaithersburg, MD, US
Contact	KAREN SIEGEL
510(k) history	7 submissions · 7 cleared · 1992-1994

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k935103/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 20, 2026