

**K935877 NIHON KOHDEN CENTRAL NURSE STATION,
CNS-8310A/CNS-8350A**

May 24, 1995
531 days to decision

K935877 · Product code: **DSI** · Cardiovascular
Source: <https://www.510kdatabase.net/k935877/>

SUBMISSION DETAILS

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Detector And Alarm, Arrhythmia (DSI)
Date received	Dec 9, 1993
Decision date	May 24, 1995
Days to decision	531 days
Third-party review	No
Summary / Statement	Summary

APPLICANT

Company	Nihon Kohden America, Inc.
Location	Foothill Ranch, CA, US
Contact	PENNI PANSELL
510(k) history	166 submissions · 163 cleared · 1979-2012

510k Database - www.510kdatabase.net

Device record: <https://www.510kdatabase.net/k935877/>. Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 24, 2026