

**K945761 LOW TEMPERATURE CAUTERY FINE TIP,  
ADJUSTABLE, REPLACEABLE TIP**Apr 10, 1995  
152 days to decisionK945761 · Product code: **HQP** · Ophthalmic  
Source: <https://www.510kdatabase.net/k945761/>**SUBMISSION DETAILS**

---

|                       |   |
|-----------------------|---|
| Decision              | Substantially Equivalent (Cleared)            |
| Submission type       | Traditional                                   |
| Device classification | Unit, Cautery, Thermal, Battery-powered (HQP) |
| Date received         | Nov 9, 1994                                   |
| Decision date         | Apr 10, 1995                                  |
| Days to decision      | 152 days                                      |
| Third-party review    | No  |
| Summary / Statement   | Statement                                     |

**APPLICANT**

---

|                |   |
|----------------|---|
| Company        | <b>Aaron Medical Industries</b>         |
| Location       | St. Petersburg, FL, US                  |
| Contact        | J R SARON                               |
| 510(k) history | 46 submissions · 46 cleared · 1990-2006 |

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k945761/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated June 29, 2026