

**K960446 MODEL 5471 & 5471L STERILE DISPOSABLE PATIENT  
CABLE ASSEMBLIES**May 22, 1996  
111 days to decisionK960446 · Product code: **DSA** · Cardiovascular  
Source: <https://www.510kdatabase.net/k960446/>**SUBMISSION DETAILS**

---

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	Cable, Transducer And Electrode, Patient, (including Connector) (DSA)
Date received	Feb 1, 1996
Decision date	May 22, 1996
Days to decision	111 days
Third-party review	No
Summary / Statement	Summary

**APPLICANT**

---

Company	<b>Medtronic Vascular</b>
Location	Walker, MI, US
Contact	VALERIE K STUHR
510(k) history	475 submissions · 453 cleared · 1977-2023

---

510k Database - [www.510kdatabase.net](http://www.510kdatabase.net)Device record: <https://www.510kdatabase.net/k960446/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).

510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 14, 2026