

**K963705 COLOR ENHANCEMENT FUNCTION FOR THE
SSA-380A**Dec 2, 1996
77 days to decisionK963705 · Product code: IYN · Radiology
Source: <https://www.510kdatabase.net/k963705/>**SUBMISSION DETAILS**

Decision	Substantially Equivalent (Cleared)
Submission type	Traditional
Device classification	System, Imaging, Pulsed Doppler, Ultrasonic (IYN)
Date received	Sep 16, 1996
Decision date	Dec 2, 1996
Days to decision	77 days
Third-party review	No
Summary / Statement	Summary

APPLICANT

Company	Toshiba America Medical Systems, In.C
Location	Tustin, CA, US
Contact	STEVEN M KAY
510(k) history	146 submissions · 146 cleared · 1989-2014

510k Database - www.510kdatabase.netDevice record: <https://www.510kdatabase.net/k963705/> Data sourced from FDA 510(k) public records (accessdata.fda.gov).510k Database is an independent platform not affiliated with the U.S. Food and Drug Administration. Always verify at accessdata.fda.gov. - Generated May 17, 2026